

# NeuroLight Lab: fNIRS Data Analysis Pipeline

Meram Mahmoud  
Department of Biomedical Engineering,  
Cairo University, Class of 2026

**Abstract**—This document describes the NeuroLight Lab project, a complete pipeline for analyzing functional Near-Infrared Spectroscopy (fNIRS) data. The pipeline includes data exploration, preprocessing, motion correction, and statistical modeling. It also provides a practical guide to handling fNIRS data with detailed examples for researchers and students.

**Index Terms**—fNIRS, Data Analysis, Motion Correction, Statistical Models, Block Averaging, GLM, LMM

## I. INTRODUCTION

Welcome to **NeuroLight Lab**, a comprehensive project that showcases a full pipeline for analyzing functional Near-Infrared Spectroscopy (fNIRS) data. This project is designed for researchers and students working with brain imaging data and provides a hands-on walkthrough from raw data to statistical modeling and interpretation.

## II. PROJECT OVERVIEW

This project demonstrates how to:

- Load and explore `.nirs` files
- Preprocess raw fNIRS signals
- Correct motion artifacts
- Apply robust statistical models, including:
  - Block Averaging
  - Generalized Linear Models (GLMs)
  - Linear Mixed Models (LMMs)

Whether you're just getting started or refining your own fNIRS analysis pipeline, this project bridges the gap between theory and practice with detailed examples.

## III. fNIRS DATA EXPLORATION

This section walks you through the process of exploring fNIRS data, including understanding the data structure, discovering the source-detector layout, calculating distances between optodes, and visualizing the raw signals.

### A. Load the Data

To begin, load the `.nirs` file using the `matfile` function in MATLAB. Once the file is loaded, check the structure to understand its contents.

The file will contain several variables, including:

- `d`: Raw light intensity (channels  $\times$  time)
- `t`: Time vector (in seconds)
- `s`: Stimulus vector (1 = event on)
- `aux`: Auxiliary signals (e.g., accelerometer data)
- `ml`: Measurement list (channel & wavelength mapping)
- `SD`: Source-detector geometry (layout and wavelengths)
- `systemInfo`: Metadata about the acquisition device


| Name       | Size     | Bytes    | Class  | Attributes |
|------------|----------|----------|--------|------------|
| SD         | 1x1      | 6773     | struct |            |
| aux        | 54675x8  | 3499200  | double |            |
| d          | 54675x28 | 12247200 | double |            |
| dStd       | 0x0      | 0        | double |            |
| ml         | 28x4     | 896      | double |            |
| s          | 54675x1  | 437400   | double |            |
| systemInfo | 1x1      | 876818   | struct |            |
| t          | 54675x1  | 437400   | double |            |
| tdml       | 28x54675 | 12247200 | double |            |

Fig. 1: file variables

For preprocessing, the main variables of interest are `d`, `t`, `s`, and `aux` for signals and events, and `SD` and `ml` for channel and optode information.


### B. Discover Optode Placement and Wavelengths

To explore the source-detector layout and wavelengths, use the following MATLAB code (see images for visualizations):



|            |     |     |
|------------|-----|-----|
| 1x2 double |     |     |
| 1          | 690 | 830 |

Fig. 2: wavelengths used



|                    |   |   |   |   |
|--------------------|---|---|---|---|
| ml x wavelengths x |   |   |   |   |
| 28x4 double        |   |   |   |   |
|                    | 1 | 2 | 3 | 4 |
| 1                  | 1 | 1 | 1 | 1 |
| 2                  | 1 | 2 | 1 | 1 |
| 3                  | 2 | 1 | 1 | 1 |
| 4                  | 2 | 2 | 1 | 1 |
| 5                  | 2 | 3 | 1 | 1 |
| 6                  | 2 | 5 | 1 | 1 |
| 7                  | 3 | 2 | 1 | 1 |
| 8                  | 3 | 4 | 1 | 1 |
| 9                  | 3 | 5 | 1 | 1 |
| 10                 | 3 | 6 | 1 | 1 |
| 11                 | 4 | 5 | 1 | 1 |
| 12                 | 4 | 6 | 1 | 1 |
| 13                 | 4 | 7 | 1 | 1 |
| 14                 | 4 | 8 | 1 | 1 |
| 15                 | 1 | 1 | 1 | 2 |
| 16                 | 1 | 2 | 1 | 2 |
| 17                 | 2 | 1 | 1 | 2 |
| 18                 | 2 | 2 | 1 | 2 |
| 19                 | 2 | 3 | 1 | 2 |
| 20                 | 2 | 5 | 1 | 2 |
| 21                 | 3 | 2 | 1 | 2 |
| 22                 | 3 | 4 | 1 | 2 |
| 23                 | 3 | 5 | 1 | 2 |
| 24                 | 3 | 6 | 1 | 2 |
| 25                 | 4 | 5 | 1 | 2 |
| 26                 | 4 | 6 | 1 | 2 |

Fig. 3: source-detector layout

### C. Get Source-Detector Distance and Visualize Optodes

Channel 1 distance: 30.02 mm  
Channel 2 distance: 30.00 mm  
Channel 3 distance: 30.00 mm  
Channel 4 distance: 30.02 mm  
...

Fig. 4: Source-detector distance.

Visualize the optodes and the layout of the fNIRS channels in 3D. This will help you understand the relative positioning of the optodes in the setup. Use visualization tools in MATLAB to achieve this.




Fig. 5: Optodes placement visualization.

## IV. PREPROCESSING STEPS

**Note:** Familiarity with digital signal processing (DSP) is recommended for understanding the preprocessing steps.

### A. Visual Inspection

Visual inspection is crucial for identifying noisy channels or time windows to exclude from analysis.




Fig. 6: subject 5

1) *Visual Inspection on subject level:* : we can discard the whole subject as the signal is corrupted.




Fig. 7: Visual inspection on subject level "sub30".

2) *Visual Inspection on channel level:* : we can discard channel due to high noise as ch2




Fig. 8: Visual inspection on channel level "sub1".

### B. Convert to Optical Density

The raw intensity data is converted to optical density using the Beer-Lambert law. This step helps to linearize the raw data and prepare it for further processing.




Fig. 9: Conversion of raw light intensity to optical density.

### C. Motion Correction

Motion artifacts are common in fNIRS data and can significantly affect signal quality. Various techniques are available to address these artifacts:

#### 1) Common Techniques for Motion Artifact Correction:

- **Wavelet Filtering:** Suppresses motion artifacts by removing outlier coefficients in the wavelet domain.
- **Principal Component Analysis (PCA):** Identifies and removes components associated with motion.
- **Spline Interpolation:** Interpolates motion-contaminated segments of the signal using smooth spline curves.
- **Kalman Filtering:** Applies a recursive filter to estimate the true signal while minimizing the impact of motion-induced noise.

#### 2) Spline-Based Motion Correction:

- Uses spline interpolation to correct segments of the signal contaminated by motion.
- Produces a corrected dataset by smoothly interpolating the affected segments and preserving the underlying signal.




Fig. 10: Spline-Based Motion Correction.

3) *Wavelet-Based Motion Correction:* Utilizes a statistical threshold to detect and eliminate motion artifacts in the wavelet-transformed data.




Fig. 11: Wavelet-based motion correction method.

### D. Filtering Techniques

Use frequency-based filters (low pass, band pass, high pass) to eliminate unwanted noise.




Fig. 12: Frequency-based filters for signal cleaning.

PCA filters can also be applied to remove high variance components.




Fig. 13: PCA filters for signal cleaning.

### E. Baseline Correction

Methods like bandpass filtering and baseline subtraction help to remove low-frequency drifts from the data.

### F. Modified Beer-Lambert Law (MBLL)

The optical density is converted to concentration changes of:

- HbO (oxyhaemoglobin)




Fig. 14: HbO

- HbR (deoxyhaemoglobin)




Fig. 15: HbR

- HbT (total haemoglobin)




Fig. 16: HbT

### G. Short Separation Channel

Use short-separation channels for GLM to regress out the component of the light source (LS) signal correlated with the short-separation (SS) signal.

## V. STATISTICAL ANALYSIS

### A. Core Questions in fNIRS Analysis

Some core questions in fNIRS analysis include:

- Which brain areas were activated?
- Are there differences between conditions or groups?

### B. Statistical Models

1) *1. Block Averaging*: Block averaging is ideal for simple experimental designs and allows for a visual inspection of hemodynamic responses before applying more complex models.

2) *2. Generalized Linear Models (GLM)*: GLMs are robust statistical models used to estimate task-related activations. The model is:

$$Y = X\beta + \epsilon$$

Where:

- $Y$  is the observed data
- $X$  is the design matrix (predictors)
- $\beta$  is the parameters to estimate
- $\epsilon$  is the error term

3) *3. Linear Mixed Models (LMMs)*: LMMs extend GLMs to account for both fixed and random effects. These models are suitable for repeated-measures or hierarchical data, and they help account for individual variability.

## VI. CONCLUSION

This pipeline provides a comprehensive approach to fNIRS data analysis, from raw data inspection to statistical modeling. By following these steps, researchers can ensure accurate analysis of hemodynamic responses and investigate brain activity reliably.

## APPENDIX

### A. fNIRS Data Preprocessing and Motion Correction

```

1 % Load subject data
2 sub_1 = matfile('S1001_run01.nirs');
3 whos(sub_1);
4
5 SD = sub_1.SD;
6 wavelengths = SD.Lambda;
7
8 ml = sub_1.ml;
9
10 % Get source-detector distance
11 SD = sub_1.SD;
12 wavelengths = SD.Lambda;
13 ml = sub_1.ml;
14 srcPos = SD.SrcPos;
15 detPos = SD.DetPos;
16 for ch = 1:4
17     s_idx = ml(ch, 1);
18     d_idx = ml(ch, 2);
19     dist = norm(srcPos(s_idx,:) - detPos(d_idx,:));
20     % Euclidean distance
21     fprintf('Channel %d distance: %.2f mm\n', ch,
22         dist);
23 end
24
25 % Visualizing the optodes
26 srcPos = SD.SrcPos;
27 detPos = SD.DetPos;
28
29 figure;
30 hold on;
31 axis equal;
32 title('fNIRS Optode Layout and Channels');
33 xlabel('X (mm)');
34 ylabel('Y (mm)');
35 zlabel('Z (mm)');
36
37 % Plot sources in red
38 for i = 1:size(srcPos, 1)
39     plot3(srcPos(i,1), srcPos(i,2), srcPos(i,3), 'ro',
40         'MarkerSize', 8, 'LineWidth', 2);

```

```

38     text(srcPos(i,1), srcPos(i,2), srcPos(i,3)+2,
          sprintf('S%d', i), 'Color', 'r');
39 end
40
41 % Plot detectors in blue
42 for i = 1:size(detPos, 1)
43     plot3(detPos(i,1), detPos(i,2), detPos(i,3), 'bo', 'MarkerSize', 8, 'LineWidth', 2);
44     text(detPos(i,1), detPos(i,2), detPos(i,3)+2,
          sprintf('D%d', i), 'Color', 'b');
45 end
46
47 % Plot channels (source-detector pairs) as lines
48 for ch = 1:size(ml, 1)
49     s_idx = ml(ch, 1);
50     d_idx = ml(ch, 2);
51
52     % Coordinates
53     s_pos = srcPos(s_idx, :);
54     d_pos = detPos(d_idx, :);
55
56     % Plot line
57     plot3([s_pos(1), d_pos(1)], [s_pos(2), d_pos(2)], [s_pos(3), d_pos(3)], 'k-', 'LineWidth', 1.5);
58
59     % Channel midpoint label
60     mid = (s_pos + d_pos) / 2;
61     text(mid(1), mid(2), mid(3), sprintf('Ch%d', ch), 'Color', 'k', 'FontSize', 8);
62 end
63
64 grid on;
65 view(3); % 3D view
66
67 % Visualize the Data
68 d = sub_1.d;
69 t = sub_1.t;
70 ml = sub_1.ml;
71 wavelengths = [690, 830]; % Known wavelengths
72
73 nChannels = size(d, 2);
74
75 figure;
76 set(gcf, 'Name', 'Raw fNIRS Data - Visual Inspection', 'NumberTitle', 'off');
77 tiledlayout('flow');
78
79 for ch = 1:nChannels
80     nexttile;
81     plot(t, d(:, ch), 'k');
82     xlabel('Time (s)');
83     ylabel('Intensity');
84
85     s = ml(ch,1);
86     d_idx = ml(ch,2);
87     lambda_idx = ml(ch,4); % Adjust if needed
88     lambda = wavelengths(lambda_idx);
89
90     title(sprintf('Ch %d (S%d-D%d, %d nm)', ch, s, d_idx, lambda));
91 end
92
93 sgtitle('Raw fNIRS Intensity Signals Across Channels');
94
95 % Plot specific channels
96 figure;
97 plot(t, d(:,1:5)); % show first 4 channels
98 xlabel('Time (s)');
99 ylabel('Intensity');
100 title('Raw Light Intensity - Channels 1 to 5');
101 legend('Ch 1', 'Ch 2', 'Ch 3', 'Ch 4', 'Ch 5');
102

```

```

103 % Convert raw intensity to optical density
104 dod = hmrIntensity2OD(d); % Homer2 function
105
106 % Plot OD for first 5 channels
107 figure;
108 plot(t, dod(:,1:5));
109 xlabel('Time (s)');
110 ylabel('Optical Density (a.u.)');
111 title('Optical Density - Channels 1 to 5');
112 legend('Ch 1', 'Ch 2', 'Ch 3', 'Ch 4', 'Ch 5');
113
114 % Wavelet-based motion correction
115 iqr_thresh = 1.5; % Interquartile range threshold ( suggested: 1.5 3 )
116
117 % Run Wavelet-based correction
118 [d_wavelet, ~] = hmrMotionCorrectWavelet(data.d, data.t, iqr_thresh);
119
120 % Replace original data with corrected data
121 data.d = d_wavelet;
122
123 % Save result
124 save('SI030_run01_WaveletCorrected.nirs', 'data');

```

Listing 1: fNIRS Data Preprocessing and Motion Correction

## B. Stimulus Visualization

```

1 % Plot the AUX channels to inspect
2 aux = finfo.aux;
3
4 figure;
5 plot(t, aux(:,1), 'r'); hold on;
6 plot(t, aux(:,2), 'b');
7 xlabel('Time (s)');
8 legend('aux1 - Stimulus onset', 'aux2 - Response cue');
9 title('AUX channels inspection');
10
11 % Detect stimulus onsets in AUX1
12 threshold = 0.5 * max(aux(:,1)); % Half the max value is safe
13 stim_onsets = find(aux(:,1) > threshold);
14
15 % Get only the rising edges
16 stim_onsets = stim_onsets([true; diff(stim_onsets) > 10]); % 10 samples apart minimum
17
18 % Convert to time
19 stim_times = t(stim_onsets);
20
21 fprintf('Detected %d stimulus onsets.\n', numel(stim_times));
22
23 % Create a new stimulus matrix "s"
24 s = zeros(length(t),1); % initialize
25
26 % Insert triggers at stimulus onset times
27 for i = 1:length(stim_onsets)
28     s(stim_onsets(i)) = 1; % mark the event
29 end
30
31 % Visualize the reconstructed stimulus
32 figure;
33 plot(t, aux(:,1)*0.5, 'r'); hold on;
34 plot(t, s, 'k');
35 legend('aux1 - Original', 's - Reconstructed stimulus');
36 xlabel('Time (s)');
37 title('Original AUX and Reconstructed Stimulus (s)');

```

Listing 2: Stimulus Detection and Visualization